

TICARCILLIN IRGASAN ANTIMICROBIC SUPPLEMENT

VIAL CONTENTS (for 250 ml of medium)

Ticarcillin.....0.25 mg
Irgasan.....0.25 mg

DIRECTIONS

Dissolve the contents of one vial with 2 ml of sterile distilled water. Prepare 250 ml of Yersinia ITC Broth Base (cat. N° 402270), autoclaved and cooled to 50°C. Add the Potassium Chlorate Supplement (cat. N° 4240065) to the basal medium together with the contents of one vial of Ticarcillin Irgasan Antimicrobial Supplement, mix well and distribute into sterile tubes.

INTENDED USE

For the enumeration of *Yersinia enterocolitica* in foods. For the details of the procedure see the technical sheet of Yersinia ITC Broth Base (cat. N° 402270).

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

Consult the material safety data sheet before the use.

PACKAGING

4240060 Ticarcillin Irgasan Antimicrobial Supplement, 10 vials, each for 250 ml of medium