

LISTERIA FRASER SUPPLEMENT (FERRIC AMMONIUM CITRATE SUPPLEMENT)

Vial contents (vial for 500 ml of medium)

Ferric Ammonium Citrate 0.25 g

Vial contents (vial for 5000 ml of medium)

Ferric Ammonium Citrate 2.5 g

DIRECTIONS

Fraser Broth- Fraser Broth Half Concentration: dissolve the contents of one vial with 5 ml of sterile distilled water under aseptic conditions. Add to 500 ml of Listeria Fraser Broth Base (cat. N° 401596) or to 500 ml of Listeria Fraser Broth Base Half Concentration (cat. N° 401594), autoclaved and cooled to 50°C. Mix well and distribute into sterile tubes or flasks.

INTENDED USE

The selective supplements described above are used for the preparation of enrichment broths and plating media for the isolation of *Listeria* spp. from food, environmental or clinical specimens. For the details of the procedure see the technical sheets of Listeria Fraser Broth Base (cat. N° 401596) or of Listeria Fraser Broth Base Half Concentration (cat. N° 401594), .

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

PACKAGING

4240056	Ferric Ammonium Citrate Supplement	10 vials, each for 500 ml of medium
42185056	Ferric Ammonium Citrate Supplement	1 vial for 5 litres of medium